

## RATIONALE:

Lung cancer screening (LCS) remains underutilized in the United States, in spite of proven mortality benefit,<sup>1,2</sup> such that in 2022 only 18% of patients who were eligible received a Low Dose CT (LDCT) scan for LCS<sup>3</sup>. We postulated that if a peripheral blood test which was highly sensitive for lung cancer were available, and more convenient for patients than a CT scan, it might augment the rate of LCS. With this aim, DELFI Diagnostics developed a blood-based genomic lung cancer screening test (pWGFRag-Lung, FirstLook Lung)<sup>4</sup> that could be obtained as a routine annual blood draw in a primary care setting. Here we present results of a cluster randomized trial to determine whether the availability of the blood test increased lung screening rates.

## METHODS

### Clinics:

Primary Care Practices across multiple US health systems were block randomized 1:1 to:

- Arm A = usual care (LDCT alone) or
- Arm B = usual care + access to the blood test

### Patients eligible for analysis:

- Met 2021 USPSTF criteria for lung cancer screening,<sup>5</sup> and
- had never completed Lung Cancer Screening (LCS) or
- were past due for their annual screening >15 months.

Outcomes were compared using time-to-event analyses.

Multi-variable analyses adjusted for patient (age, sex, race/ethnicity, pack-year smoking history, current smoking status) and medical practice (network, rural/urban status) characteristics.

Patients were followed up to 11 months.

### Outcomes:

1. **Primary Outcome** = completion of "any lung cancer screening" (LDCT or blood test)
2. **Secondary Outcome** = completion of LDCT.

## RESULTS

28 primary care practices were randomized:

15 practices to usual care (Group A; 1662 eligible patients), and

13 practices to usual care plus access to the blood test (Group B; 1291 eligible patients).

**Primary outcome:** Any lung cancer screening (LDCT or blood test) occurred more frequently in Group B practices (20.9%) than Group A practices (8.4%): multi-variable HR = 2.73 (95% CI: 1.91, 3.89; p < 0.001).


**Secondary outcome:** LDCT screening was more common in Group B (13.4%) than Group A (8.4%): multi-variable HR: 1.47 (95% CI: 1.14, 1.89; p < 0.01).

## REFERENCES:

- 1: "Reduced Lung-Cancer Mortality with Low-Dose Computed Tomographic Screening": The NLST Research Team. N Engl J Med 2011;365:395-409, VOL. 365 NO. 5
- 2: "Reduced Lung-Cancer Mortality with Volume CT Screening in a Randomized Trial": Harry J. de Koning, M.D et al., N Engl J Med 2020;382:503-513, VOL. 382 NO. 6
- 3: The American Lung Association's "State of Lung Cancer 2025"
- 4: "Clinical Validation of a Cell-Free DNA Fragmentome Assay for Augmentation of Lung Cancer Early Detection", Mazzone et al Cancer Discov. 2024 Jun 6;14(11):2224-2242.
- 5: USPSTF. JAMA. 2021;325(10):962-970.
- 6: Presented at ATS 2026; Poster Board #608: "FIRSTLUNG L301: Baseline Socioeconomic Characteristics of a Cluster Randomized Clinical Utility Trial Evaluating Blood-Based Lung Cancer" Screening Across Diverse Communities.
- 7: Kazi MAI, Qureshi I, Sharma A, et al. From detection to delay: real-world gaps in post-Colorectal colonoscopy adherence. J Gastrointest Cancer. 2025;56(1). doi:10.1007/s12029-025-01350-5

**\*DISCLAIMER**  
 The FirstLook Lung test is a laboratory-developed test. This test was developed, and its performance characteristics were determined by DELFI Diagnostics. It has not been cleared or approved by the US Food and Drug Administration (FDA). The laboratory is regulated under the Clinical Laboratory Improvement Act (CLIA) as qualified to perform high-complexity clinical tests.

**Figure 1: USPSTF Recommend Cancer Screening and Association with Uptake**




**Figure 2: Performance of the Blood-based Test made Available in the Study = pWGFRag-Lung (FirstLook Lung)<sup>4</sup>**

Screening-Population	Performance
Screening-Population Sensitivity <sup>1</sup>	80%
Screening-Population Specificity <sup>2</sup>	56%
NPV <sup>1</sup>	99.8%

<sup>1</sup>Test performance adjusted to the USPSTF eligible lung cancer screening population.  
<sup>2</sup>Unpublished data on file DELFI Diagnostics Inc.

**Figure 3: Cluster Randomized Trial design of FIRSTLUNG**  
 CME = Continuing Medical Education



**Table 1: Clinic Characteristics Balanced Between Arms**

Geographic regions are largely urban; enrollment clinics were selected for catchment areas of high socioeconomic diversity<sup>6</sup>

Clinic and Patient Characteristics Overall and by Analytic Allocation	Arm Assignment		
	Overall, N = 28	Usual Care (Arm A) <sup>6</sup> , N = 15	Intervention (Arm B), N = 13
Health Network, n (%)			
Network 1	20 (71%)	10 (67%)	10 (77%)
Network 2	8 (29%)	5 (33%)	3 (23%)
Geographic Region, n (%)			
Urban	26 (93%)	13 (87%)	13 (100%)
Rural	2 (7%)	2 (13%)	0 (0%)
Patients per Clinic, Median (IQR)	85.5 (52.5 - 141.8)	89.0 (64.5 - 133.5)	79.0 (34.0 - 139.0)
Clinicians per Clinic, Median (IQR)	8.0 (4.8 - 11.0)	7.0 (4.5 - 10.5)	9.0 (6.0 - 11.0)

**Table 2: Patient Lung Cancer Screening History Overall and by Analytic Allocation: Balanced Between Arms**

The Study Included Only Patients That Were Behind On Screening

Patient Characteristics Overall and by Analytic Allocation	Arm Assignment		
	Overall, N = 2,953	Usual Care (Arm A) <sup>6</sup> , N = 1,662	Intervention (Arm B), N = 1,291
LDCT Screening History			
Behind on Screening (>15 months)	222 (7.5%)	134 (8.1%)	88 (6.8%)
Never Screened	2,414 (81.7%)	1,369 (82.4%)	1,045 (80.9%)
Newly Eligible, Age criteria	59 (2.0%)	31 (1.9%)	28 (2.2%)
Unknown/Missing	258 (8.7%)	128 (7.7%)	130 (10.1%)


**Table 3: Patient Characteristics Overall and by Analytic Allocation. Enrollment of Clinical Sites achieved a diverse screening population**

Q1 = Interquartile range SD = Standard deviation


Patient Characteristics Overall and by Analytic Allocation	Arm Assignment		
	Overall, N = 2,953	Usual Care (Arm A) <sup>6</sup> , N = 1,662	Intervention (Arm B), N = 1,291
Age, years, Mean (SD)	64.5 (7.8)	64.6 (7.8)	64.3 (7.9)
Sex, n (%)			
Female	1,512 (51.2%)	844 (50.8%)	668 (51.7%)
Male	1,441 (48.8%)	818 (49.2%)	623 (48.3%)
Race/Ethnicity <sup>6</sup> , n (%)			
Non-Hispanic White	2,023 (68.5%)	1,188 (71.5%)	835 (64.7%)
Non-Hispanic Black	708 (24.0%)	351 (21.1%)	357 (27.7%)
Hispanic	96 (3.3%)	62 (3.7%)	34 (2.6%)
Other	126 (4.3%)	61 (3.7%)	65 (5.0%)
Smoking Status, n (%)			
Current	1,500 (50.8%)	836 (50.3%)	664 (51.4%)
Former	1,453 (49.2%)	826 (49.7%)	627 (48.6%)
Pack-Years, Median, (IQR)	34.0 (25.0 - 45.6)	35.0 (25.0 - 46.2)	33.8 (25.0 - 45.0)
Years Since Quit (Former smokers), Mean (SD)	7.5 (4.3)	7.5 (4.3)	7.5 (4.3)
Chronic Obstructive Pulmonary Disease (COPD), n (%)	1,368 (46.3%)	770 (46.3%)	598 (46.3%)
Cardiovascular Disease, n (%)	1,361 (46.1%)	760 (45.7%)	601 (46.6%)
Asthma, n (%)	567 (19.2%)	295 (17.7%)	272 (21.1%)

**Abbreviations**  
 CI = Confidence Interval  
 CME = Continuing Medical Education  
 COPD = Chronic Obstructive Pulmonary Disease  
 HR = Hazard Ratio  
 IQR = Interquartile range  
 LCS = Lung Cancer Screening  
 LDCT = Low-dose computed tomography  
 NH = Non-Hispanic  
 NPV = negative predictive value  
 SD = Standard deviation  
 USPSTF = United States Preventive Services Task Force

**Figure 6: Blood-Based Testing Resulted in Clinically Meaningful Increase in Any Lung Cancer Screening (LCS) and Screening by LDCT**



**Figure 7: Observed and Adjusted Effect Sizes Across Subgroups of Interest**



Abbreviations: LCS, lung cancer screening; CI, control arm; HR, hazard ratio; CI, confidence interval; HR, hazard ratio; CME, continuing medical education; RMST, restricted mean survival time; n, size; LDCT, low-dose CT; NH, Non-Hispanic; COPD, Chronic Obstructive Pulmonary Disease.

## Availability of Blood-based Screening increased Lung Cancer Screening rates with a modeled 5 year potential to avert more deaths

Primary Outcome: Any Lung Cancer Screening

**2.73X**

Screening Rate Increase

Any LCS: 20.9% (intervention) vs. 8.4% (usual care)

Adjusted HR 2.73; P < .001

Secondary Outcome: LDCT Screening

**1.47X**

LDCT Rate Increase

LDCT: 13.4% (intervention) vs. 8.4% (usual care)

Adjusted HR 1.47; P < .003

LDCT Follow-up after Blood-based Testing

Post-positive result LDCT follow-up

**42.4%**

Post-negative result LDCT follow-up

**5.6%**

Projected Population Health Impact

**~3,178**


Lung Cancer Deaths Preventable

Projected over 5 years, among 13 million simulated individuals, with 3% annual uptake increase (1,896-4,461)

Potential to avert an average of ~3,178 more deaths than LDCT only\*

\*28 clinics and 2,953 patients who were behind on lung screening recommendations  
 \*LDCT only: averted average of 5,980 deaths compared to No Screening  
 Any LCS: averted average of 9,158 deaths compared to No Screening

**Figure 8: Cumulative Incidence of LDCT Following Return of a Positive Blood-Test Result, Stratified by Timing of Blood Test Result**



Data indicates that there is a "learning curve" with higher follow through in the 2nd half of the study.

Estimated cumulative incidence: Early Adopter (blue line), 38.2% (95% CI, 19.5% - 52.6%); Late Adopter (red line), 48.2% (95% CI, 20.8% - 66.1%). 95% CIs derived via cluster bootstrap resampling.


Overall cumulative incidence of LDCT following return of a positive blood-based screening test, stratified by Result:  
 Negative blood test → 5.6% (95% CI, 1.5% - 13.0%) go on to LDCT  
 Positive blood test → 42.4% (95% CI, 31.2% - 55.9%) go on to LDCT

These rates are similar to the early data on colonoscopy rates following positive stool-based testing results for colorectal cancer screening.<sup>8</sup>

**Figure 9: Projected 5 year Population Health Impact of LCS Blood-based Test Availability in Primary Care Clinics**

Impact of blood-based screening vs LDCT only:

- Increases number of cancers detected
- Increases percent of Stage I detected cancers
- Decreases percent of stage IV lung cancers
- Averts ~3,178 additional deaths due to lung cancer



Modeled with 0.7% lung cancer prevalence, 24% LDCT false-positive rate,<sup>7</sup> 80% and 56% LCS blood test sensitivity and specificity. 3.0% hypothetical annual increase for both usual care and blood-based cancer screening. Method: The model is previously described in Mazzone et al<sup>8</sup>. Trial results were applied in Monte Carlo simulations to a base population of the 13 million currently eligible but unscreened or underscreened Americans<sup>3</sup>. Outcomes, including lung cancer screenings, lung cancer diagnosis by stage, and deaths due to lung cancer, were probabilistically generated. The population was cycled and aged in annual increments over 5 years.

## CONCLUSIONS

In this Cluster randomized trial of 28 clinics and 2,953 patients who were behind on lung screening recommendations, blood-based screening meaningfully increased Lung Cancer Screening rates. A 5 year modeled outcomes analysis reveals potential to avert ~3178 more deaths than by LDCT only.

1. The availability of a blood-based test significantly increased overall screening rates (20.9%) and low-dose CT rates (13.4%) compared with usual care (8.4%).
2. The effect is consistent across subgroup stratifications.
3. Post-positive clinical follow-up of 42.4% without dedicated navigation programs, sustained education or a focus on program efficiencies.
  - Similar to rates of follow-up during early implementation of stool-based testing for Colorectal Cancer screening.<sup>8</sup>
4. Conservative 5 year outcomes modeling suggests that an average of 3,178 (1,896 - 4,461) more deaths due to lung cancer can be averted with access to this blood-based testing approach compared to LDCT only, assuming stable or modestly increasing screening uptake, respectively.

Given the low rates of lung cancer screening via LDCT, access to a lung cancer screening blood test, such as the pWGFRag-Lung test, could substantially increase LCS rates and reduce lung cancer mortality.